# **TECHNICAL ADVANCE**

**Open Access** 

# Preliminary experience of stent-assisted coiling of wide-necked intracranial aneurysms with a single microcatheter



Keun Young Park<sup>1</sup>, Chang Ki Jang<sup>1</sup>, Jae Whan Lee<sup>1</sup>, Dong Joon Kim<sup>2</sup>, Byung Moon Kim<sup>2</sup> and Joonho Chung<sup>1,3\*</sup>

# **Abstract**

**Background:** The purpose of this study was to report our preliminary experience of stent-assisted coiling (SAC) of wide-necked intracranial aneurysms with a single microcatheter in patients with parent arteries that were small-caliber, with stenosis, or a very tortuous course.

**Methods:** Between March 2018 and December 2018, we treated 394 aneurysms in 359 patients with endovascular treatment. Among 197 aneurysms treated by SAC, there were 16 cases (all wide-necked unruptured aneurysms) treated by SAC with a single microcatheter and a Neuroform Atlas stent. Follow-up angiography was performed at 6 to 12 months after SAC, and clinical follow-up was performed from 6 to 12 months in all patients.

**Results:** The reasons for SAC with a single 0.0165-in. microcatheter were small-caliber (n = 4), stenosis (n = 2), and very tortuous course (n = 10) of the parent arteries. There was no complication related to delivering or deploying the Neuroform Atlas stent as well as no failure of selecting aneurysm by cell-through technique. All patients had a modified Rankin score of 0 at discharge and at follow-up. Initial angiographic results showed six cases (37.5%) of complete occlusion. In follow-up angiographies, 12 cases (75.0%) achieved compete occlusion.

**Conclusion:** When performing SAC of wide-necked intracranial aneurysms in parent arteries with small-caliber, stenosis, or a very tortuous course, cell-through SAC using a single microcatheter and a Neuroform Atlas stent within a 5 Fr- (or smaller) guiding or intermediate catheter might be a useful option.

Keywords: Endovascular treatment, Intracranial aneurysm, Neuroform atlas, Stent assisted coiling, Stent through

# **Background**

Despite accumulated experience and improved understanding of the devices, endovascular treatment of widenecked intracranial aneurysms continues to be difficult to perform due to limitations arising from cerebrovascular anatomical structures [1]. When performing stent-assisted coiling (SAC), a 6-Fr or larger guiding catheter is usually needed to support two microcatheters; one for stent delivery and the other for coil delivery. A guiding catheter is placed in the internal carotid artery (ICA) for anterior circulation aneurysms and in the vertebral artery (VA) for posterior circulation aneurysms. However,

when placing a 6-Fr guiding catheter it is possible to encounter ICAs or VAs that are small-caliber, exhibit stenosis, or a very tortuous course. In such cases, a smaller (5-Fr or less) guiding catheter or an intermediate catheter can be placed to support a microcatheter, but not two microcatheters. Thus, SAC for wide-necked intracranial aneurysms cannot be performed with two microcatheters if they have proximal parent arteries exhibiting a small-caliber, stenosis, or a very tortuous course.

The purpose of this study was to report our preliminary experience with cell-through SAC for wide-necked intracranial aneurysms with a single 0.0165-in. microcatheter in patients with parent arteries of small-caliber, with stenosis, or exhibiting a very tortuous course.

Full list of author information is available at the end of the article



<sup>\*</sup> Correspondence: ns.joonho.chung@gmail.com

<sup>&</sup>lt;sup>1</sup>Department of Neurosurgery, Severance Hospital, Yonsei University College of Medicine, 50-1, Yonsei-ro, Seodaemun-gu, Seoul 03722, Republic of Korea <sup>3</sup>Severance Institute for Vascular and Metabolic Research, Yonsei University College of Medicine, Seoul, Republic of Korea

Park et al. BMC Neurology (2019) 19:245 Page 2 of 9

#### Methods

# Patient selection

This retrospective study was approved by our Institutional Review Board, and the requirement for informed consent was waived. From March 2018, the Neuroform Atlas (Stryker Neurovascular, Fremont, USA) has been approved and used in our country. Since then, we treated 394 aneurysms (336 unruptured and 58 ruptured) in 359 patients with endovascular treatment. Among 197 aneurysms treated with SAC, there were 16 cases (all wide-necked unruptured aneurysms) treated with cell-through SAC using a single microcatheter, Excelsior SL-10 (Stryker Neurovascular). Among the 58 ruptured aneurysms during the same period, we tried to perform endovascular treatment other than using stents due to our results of stent-assisted coiling in ruptured aneurysms [2].

All the 16 cases underwent digital subtraction angiography (DSA) and rotational angiography (Philips Allura FD20 Clarity System, Philips Medical Systems, Best, The Netherlands) with three-dimensional reconstructions to characterize aneurysm and parent artery anatomy (Allura 3D-RA workstation, Philips Medical Systems). From these DSA results, proximal parent arteries, including the extracranial ICA or VA, were determined to use a smaller (< 6-Fr) guiding catheter due to its small-caliber, stenosis, or tortuous course. Wide-necked aneurysms were defined as having a neck width  $\geq 4\,\mathrm{mm}$  or a dome-to-neck ratio of < 2 [2, 3].

# **Endovascular treatment**

The treatment decision was made by agreement of microvascular neurosurgeons and neurointerventionalists. All SAC were performed under general anesthesia. A 6-Fr femoral sheath was inserted into the right common femoral artery. A 5-Fr guiding catheter Envoy (Codman Neurovascular, Miami Lakes, FL, USA) was positioned in the ICA or VA. In some cases, a 6-Fr Shuttle (Cook Medical, Bloomington, IN, USA) was positioned in the common carotid artery or the subclavian artery followed by a 5-Fr intermediate catheter (Sofia, MicroVention, Tustin, CA, USA) in the ICA or VA. Intravenous systemic heparin (50 U/kg) was administered while placing the guiding catheter. Patients took a daily dose of 75 mg of clopidogrel and 100 mg of aspirin for more than 7 days before the SAC procedure.

An Excelsior SL-10 microcatheter was navigated to parent vessels to deliver the Neuroform Atlas stent to cover the wide-neck of the aneurysm (Fig. 1a). After deploying the Neuroform Atlas (Fig. 1b), the microcatheter was placed just proximal to the Neuroform Atlas (Fig. 1c). After removing the pusher wire of the Neuroform Atlas, a microwire was advanced through the stent proximal closed-cell segment. When the microwire

could not be advanced through the stent or was stuck by the struts, a 0.014-in. microwire with its tip flexed into a loop (microwire looping technique) provided easy navigation through the stent (Fig. 1d) and was placed just distal to the aneurysm neck (Fig. 1e). Using torque, we advanced the tip of the microwire point to the aneurysm neck and pulled the microwire back slowly. Then, the tip was caught by the cell of the stent and was inserted into the aneurysm sac (Fig. 1f). Following the microwire, the microcatheter was advanced to the neck of the aneurysm (Fig. 1g). When the microcatheter could not be advanced through the stent or was stuck in the struts, pushing the microcatheter very gently against the strut while torqueing and withdrawing the microwire simultaneously allowed the microcatheter to overcome the thin-flexible struts and advance into the aneurysm (Fig. 1h). Then, coil embolization with the cell-through technique was conducted.

After the procedure, patients were prescribed 75 mg of clopidogrel daily for 3 months and 100 mg of aspirin daily for at least 12 months. Magnetic resonance imaging was routinely performed within 24 h after SAC for all patients and included diffusion-weighted imaging, T2-weighted imaging, fluid-attenuated inversion recovery, and gradient echo imaging.


# Clinical and radiographic outcomes

Patient and aneurysm characteristics, procedure-related complications, angiographic results, and clinical outcomes were reviewed retrospectively. Procedure-related events were defined as any (symptomatic or asymptomatic) events during the procedure. Immediate postcoiling and follow-up angiographic results were evaluated by 2 independent investigators who were not involved in patient care and categorized according to the modified Raymond-Roy classification (Class I: complete obliteration, Class II: residual neck, Class IIIa: residual aneurysm: contrast within coil interstices, Class IIIb: residual aneurysm: contrast along aneurysm wall) [4]. Clinical outcomes were assessed using the modified Rankin Scale (mRS) and evaluated at the time of discharge and out-patient follow-up by an independent investigator who was not involved in patient care.

# Statistical analysis

All statistical analyses were performed using IBM SPSS Statistics version 22.0 (IBM, Armonk, New York, USA) in consultation with a biostatistician. Mann-Whitney U test was used for numeric variables. Chi-square test or Fisher exact test were used for nominal variables. A *P* value less than 0.05 for a 95% confidence interval was considered statistically significant.

Park et al. BMC Neurology (2019) 19:245 Page 3 of 9



**Fig. 1 a** A microcatheter, Excelsior SL-10, was navigated to parent vessels for delivering the Neuroform Atlas stent to cover the wide-neck of the aneurysm. **b** The Neuroform Atlas stent was being deployed. **c** The microcatheter was placed just proximal to the Neuroform Atlas deployed. **d** When the microwire could not be advanced through the stent or was stuck by the struts, a 0.014-in. microwire using its tip flexed into a loop (microwire looping technique) provided easy navigation through the stent, and (**e**) placed just distal to the aneurysm neck. **f** Make the tip of the microwire point to the aneurysm neck by using torque and pull the microwire back slowly. Then the tip can be caught by the cell of the stent and inserted into the aneurysm sac. **g** Following the microwire, the microcatheter was advanced to the neck of the aneurysm. **h** When the microcatheter could be advanced through the stent or was stuck in the struts, the microcatheter could overcome the thin-flexible struts and be advanced into the aneurysm by pushing the microcatheter very gently against the strut while torqueing and withdrawing the microwire simultaneously

# **Results**

The data and outcomes from all 16 cases are summarized in Table 1. There were 4 men and 12 women. The ages of the patients ranged from 50 to 74 years, with a mean of 62.8 years. Six aneurysms were located in the anterior communicating artery, five were in the middle cerebral artery, and five were at the basilar bifurcation. All 16 aneurysms were classified as small (< 10 mm) with a mean maximal length of 4.9 mm and a mean neck size of 4.0 mm. The reasons for using the cell-through technique were parent artery small-caliber (n = 4), stenosis (n = 2), or tortuous course (n = 10). There were no complications related to delivering or deploying the Neuroform Atlas stent as well as no failures in selecting the aneurysm by cell-through technique. However, two asymptomatic procedure-related events, left proximal VA dissections, occurred during the placement of a guiding catheter (Case 5 & 10).

Magnetic resonance angiography (MRA) or DSA were performed in all patients at 6 to 12 months after SAC, and clinical follow-up was performed from 6 to 12

months (mean 8.7 months) after SAC. All patients had mRS 0 at discharge and at follow-up. Initial angiographic results according to the Raymond-Roy classification showed six cases (37.5%) of class I, seven cases (43.8%) of class II, and three cases (18.8%) of class IIIa. In follow-up angiographies, 12 cases (75.0%) showed class I, three cases (18.8%) showed class II, and one case (6.3%) showed class IIIa.

Comparing to 181 cases undergone SAC using standard two microcatheters during the same period, there were no significant differences in baseline characteristics, procedure-related complications (asymptomatic and symptomatic), and angiographic results (Table 2).

# Case illustration

#### Case 5

A 55-year-old woman presented with a long-standing history of headache. Her computed tomography (CT) and CT angiography revealed no intracranial hemorrhage and a basilar bifurcation aneurysm. DSA demonstrated a widenecked (4.1 mm) basilar bifurcation aneurysm sized 4.9

Park et al. BMC Neurology (2019) 19:245 Page 4 of 9

**Table 1** All the data and outcomes of patients

Pt.	Sex	Age	Aneurysm location	Length (mm)	Neck size (mm)	Reason for cell-through technique	Stent size (mm)	Procedure- related event	Initial angiographic results	mRS at discharge	Follow-up period (months)	Follow-up angiographic results	mRS at follow-up
1	F	62	MCA bif., Left	4.3	3.4	Tortuosity	3.0 × 15	No	1	0	12	I	0
2	F	70	Acom, left	4.4	3.3	Stenosis	$3.0 \times 15$	No	II	0	12	1	0
3	F	69	M1, left	3.9	3.8	Tortuosity	$4.0 \times 21$	No	Illa	0	11	II	0
4	F	67	Acom, right	5.1	3.9	Tortuosity	3.0 × 15	No	1	0	11	I	0
5	F	55	BA bif.	4.9	4.1	Small-caliber	$4.0 \times 21$	Yes	II	0	10	1	0
6	Μ	50	Acom, left	4.9	3.7	Small-caliber	$3.0 \times 15$	No	Illa	0	10	Illa	0
7	F	60	M1, left	4.6	3.9	Stenosis	$3.0 \times 15$	No	II	0	10	1	0
8	М	64	Acom, right	4.8	3.9	Small-caliber	3.0 × 15	No	1	0	9	I	0
9	F	58	BA bif.	5.2	4.1	Tortuosity	$4.0 \times 21$	No	II	0	8	1	0
10	F	61	BA bif.	4.4	3.5	Tortuosity	$4.0 \times 21$	Yes	1	0	8	1	0
11	F	58	BA bif.	5.8	4.0	Tortuosity	$4.0 \times 21$	No	II	0	8	1	0
12	F	69	BA bif,	8.3	5.5	Small-caliber	$4.0 \times 21$	No	1	0	6	1	0
13	F	74	Acom, left	4.7	3.0	Tortuosity	$3.0 \times 15$	No	1	0	6	1	0
14	Μ	64	A1, left	4.2	4.7	Tortuosity	$3.0 \times 15$	No	II	0	6	II	0
15	F	58	M1, left	4.4	4.1	Tortuosity	$3.0 \times 15$	No	II	0	6	1	0
16	Μ	65	M1, right	4.6	4.8	Tortuosity	$3.0 \times 15$	No	Illa	0	6	II	0

Acom Anterior communicating artery, BA Basilar artery, bif. Bifurcation, CO Complete occlusion, F Female, M Male, MCA Middle cerebral artery, mRS Modified Rankin Scale, NA Not applicable, NR Neck remnant, PO Partial occlusion Angiographic results were classified according to the Raymond-Roy classification

mm with a daughter sac on left lateral side of the sac (Fig. 2a). On DSA, left proximal VA was dominant compared to right VA (Fig. 2b and c) so that we decided to select left VA for a guiding catheter during endovascular treatment. However, a 6-Fr guiding catheter could not pass through a left VA tortuosity at the proximal segment (a white arrowhead in Fig. 2d). Unfortunately, there was asymptomatic dissection of the left VA (a white arrow in Fig. 2d). We changed our plan for using a 5Frguiding catheter in the left VA and selected the right VA, which has a smaller diameter compared to the left VA. A single microcatheter (Excelsior SL-10 straight) was navigated to the right posterior cerebral artery and we deployed a Neuroform Atlas stent (4.0 mm × 21 mm) from the right posterior cerebral artery to the basilar artery (black arrowheads in Fig. 2e). After deploying the stent, the microcatheter was placed just proximal to the deployed stent (a black arrow in Fig. 2e). After removing the stent pusher wire, a 0.014-microwire and the same microcatheter were navigated easily through the stent to the aneurysm (Fig. 2f). Coiling was performed (Fig. 2g) and initial angiographic results showed a neck remnant of the aneurysm (Fig. 2h). On the 6-month follow-up angiography, the aneurysm was completely occluded (Fig. 2i).

#### Case 7

A 60-year-old woman experienced headaches, and a left middle cerebral artery aneurysm was diagnosed on MRA. She was willing to treat the lesion aggressively because of a family history of ruptured intracranial aneurysms. DSA revealed a middle cerebral artery aneurysm approximately 4.6 mm in maximal length with a neck size of 3.9 mm and carotid stenosis of about 40% in the left carotid bulb. For endovascular treatment of the aneurysm, a 5Fr-guiding catheter was placed in the left ICA, but it seemed to be stuck in the stenotic segment of the ICA (white dots inside of a white circle in Fig. 3a). A single microcatheter (pre-shaped Excelsior SL-10 45 degree) was steam-shaped into a "Z" shape and navigated to the left MCA. We deployed a Neuroform Atlas stent  $(3.0 \text{ mm} \times 15 \text{ mm})$  in left M1 (black arrowheads) and selected the aneurysm with the same microcatheter using the cell-through technique (Fig. 3b). A neck remnant resulted in the final angiography after coiling (Fig. 3c). On 6-month follow-up angiography, the aneurysm was completely occluded (Fig. 3d).

# Case 9

A 58-year-old woman with sudden bursting severe headache was admitted. She was alert and there was no intracranial hemorrhage on her brain CT and MRI. DSA Park et al. BMC Neurology (2019) 19:245 Page 5 of 9

**Table 2** Comparison of baseline characteristics, procedure-related complications, and angiographic outcomes between 16 cases undergone Neuroform Atlas stent-assisted coiling using a single microcatheter and 181 cases undergone stent-assisted coiling using two or more microcatheters

Variables	Atlas using a single microcatheter ( $n = 16$ )	SAC using two or more microcatheters (n = 181)	P value	
Age (mean ± SD, years)	62.8 ± 6.2	59.4 ± 10.2	0.525ª	
Female (n, %)	12 (75.0)	136 (75.1)	0.913	
Hypertension (n, %)	6 (37.5)	69 (38.1)	0.846	
Diabetes (n, %)	2 (12.5)	21 (11.6)	0.782	
Smoking (n, %)	2 (12.5)	26 (14.4)	0.636	
Aneurysm location (n, %)			0.253	
Anterior circulation	11 (68.8)	145 (80.1)		
Posterior circulation	5 (31.3)	36 (19.9)		
Neck size (mean ± SD, mm)	$4.0 \pm 0.6$	$4.5 \pm 2.5$	0.225 <sup>a</sup>	
Aneurysm size (mean $\pm$ SD, mm)	$4.9 \pm 1.0$	$5.9 \pm 2.8$	0.175 <sup>a</sup>	
Stent used (n, %)			0.061	
Neuroform Atlas	16 (100)	120 (66.3)		
Enterprise2	0 (0)	26 (14.4)		
LIVS or LVIS Jr.	0 (0)	32 (17.7)		
Procedure-related complications (n, %)			0.112	
Asymptomatic	2 (12.5)	9 (5.0)		
Symptomatic	0	1 (0.6)		
Initial angiographic results (n, %)			0.088	
Raymond-Roy I	6 (37.5)	82 (45.3)		
Raymond-Roy II	7 (43.8)	64 (35.4)		
Raymond-Roy III	3 (18.8)	35 (19.3)		
Follow-up angiographic results (n, %)	n = 16	n = 155	0.814	
Raymond-Roy I	12 (75.0)	115 (74.2)		
Raymond-Roy II	3 (18.8)	31 (20.0)		
Raymond-Roy III	1 (6.3)	9 (5.8)		

<sup>&</sup>lt;sup>a</sup>Mann-Whitney U test. SAC Stent-assisted coiling, SD Standard deviation


showed a wide-necked (4.1 mm) basilar artery bifurcation aneurysm sized 5.2 mm with a speculate in the left posterior-lateral direction (Fig. 4a). We planned to treat the aneurysm by using SAC. On DSA, the left proximal VA was dominant compared to the right VA (Fig. 4b and c) so we decided to select the left VA for a guiding catheter during endovascular treatment. Neither a 6-Fr nor 5-Fr guiding catheter went up through the left proximal VA (a black arrow) tortuous segment. We placed a 5-Fr guiding catheter in the left proximal VA and a single microcatheter (preshaped Excelsior SL-10 45 degree) was navigated from the VA to the right posterior cerebral artery (Fig. 4d). We deployed a Neuroform Atlas stent  $(4.0 \text{ mm} \times 21)$ mm) from the right posterior cerebral artery to the basilar artery (black arrowheads). We selected the aneurysm with the same microcatheter through the stent struts and coiling was performed (Fig. 4e).

Initial angiographic results showed a neck remnant of the aneurysm (Fig. 4f).

# **Discussion**

In the present study, we described 16 cases of SAC treated with the cell-through technique using a single microcatheter (Excelsior SL-10) and a Neuroform Atlas stent in patients with unruptured wide-necked intracranial aneurysms. There was no complication related to delivering or deploying the Neuroform Atlas stent as well as no failure of selecting aneurysms by the cell-through technique. Not the cell-through technique but the jailing technique is our routine method for SAC. However, we found that this cell-through technique might be feasible in patients with small-caliber, stenotic, or very tortuous parent arteries.


Recently, complicated wide-necked aneurysms rather than simple saccular aneurysms have been treated by Park et al. BMC Neurology (2019) 19:245 Page 6 of 9



**Fig. 2 a** 3D-DSA demonstrated a wide-necked (4.1 mm) basilar bifurcation aneurysm sized 4.9 mm with a daughter sac on the left lateral side of the sac. **b** and **c** Left proximal vertebral artery (VA) was dominant compared to the right VA so that we decided to select the left VA for a guiding catheter during endovascular treatment. **d** A 6-Fr guiding catheter could not pass through the tortuosity of the left VA at the proximal segment (a white arrowhead). There was asymptomatic dissection of the left VA (a white arrow). **e** With a 5Fr-guinding catheter in the right VA, a single microcatheter (Excelsior SL-10 straight) was navigated to the right posterior cerebral artery and we deployed a Neuroform Atlas stent (3.0 mm × 21 mm) from the right posterior cerebral artery to the basilar artery (black arrowheads). After deploying the stent, the microcatheter was placed just proximal to the stent deployed (a black arrow). **f** A 0.014-microwire and the same microcatheter were navigated easily through the stent to the aneurysm. **g** Coiling was performed. **h** The initial angiographic result showed a neck remnant of the aneurysm. (i) On 6-month follow-up angiography, the aneurysm was completely occluded


endovascular methods, but surgeons need more than a single microcatheter to treat those aneurysms. Multiple microcatheter, balloon-assisted, or stent-assisted techniques are needed. For that purpose, surgeons also need a strong and firm guiding catheter to support those devices. A guiding catheter has a very important role in supporting devices being delivered to tortuous intracranial vessels when performing endovascular treatments, including microcatheters, coils, balloons, or stents [5–7]. However, when faced with an ICA or VA of small-

caliber, with stenosis, or a very tortuous course, it can be difficult to place a 6Fr- or larger guiding catheter to perform multiple microcatheter, balloon-assisted, or stent-assisted techniques. In such cases, a 5Fr- (or smaller) guiding catheter or intermediate catheter would be placed and there will be only one option to use a single microcatheter. Thus, in those situations, we tried to perform SAC with a cell-through technique for widenecked intracranial aneurysms by using a single microcatheter and a Neuroform Atlas stent.



**Fig. 3 a** A 5Fr-guiding catheter was placed in the left internal carotid artery (ICA) with stenosis of about 40% in the left carotid bulb. The guiding catheter seemed to be stuck in the stenotic segment of the ICA (white dots inside of a white circle). **b** A single microcatheter (pre-shaped Excelsior SL-10 45 degree) was steam-shaped into a "Z" shape and navigated to the left MCA. We deployed a Neuroform Atlas stent (3.0 mm × 15 mm) in the left M1 (black arrowheads) and selected the aneurysm with the same microcatheter using the cell-through technique. **c** A neck remnant was seen in the final angiography after coiling. (**d**) On 6-month follow-up angiography, the aneurysm was completely occluded

Park et al. BMC Neurology (2019) 19:245 Page 7 of 9



**Fig. 4 a** 3D-DSA showed a basilar artery bifurcation aneurysm sized 5.2 mm with a spiculate on the left posterior-lateral direction. **b** and **c** Left proximal vertebral artery (VA) was dominant compared to the right VA so that we decided to select the left VA for a guiding catheter. Neither a 6-Fr nor 5-Fr guiding catheter went up through the tortuous segment of the left proximal VA (a black arrow). **d** We placed a 5-Fr guiding catheter in the left proximal VA and a single microcatheter (pre-shaped Excelsior SL-10 45 degree) was navigated from the VA to the right posterior cerebral artery. **e** We deployed a Neuroform Atlas stent (3.0 mm × 21 mm) from the right posterior cerebral artery to the basilar artery (black arrowheads). We selected the aneurysm with the same microcatheter through the stent struts and coiling was performed. **f** Initial angiographic results showed a neck remnant of the aneurysm

The Neuroform Atlas stent is the latest version of the Neuroform series reinforcing its navigability in a smaller microcatheter, stability to the vessel wall (high radial force), accurate placement (ease of use with very low foreshortening), wall apposition, conformability, and coil protection. The stent can be delivered via a 0.0165-in. microcatheter, the smallest profile, allowing surgeons to perform SAC with a single microcatheter, a sequential technique where the same microcatheter that is used to place the stent is then used to coil the aneurysm. The cell size of the Neuroform Atlas stent has been reduced compared to that of its predecessor, the Neuroform stent, to achieve better coil protection in aneurysms and to allow the use of smaller coils. Additionally, the width of struts has been reduced so that the flexibility and

conformability of the stent are improved. Furthermore, the segmental opening of the struts allows stable positioning and direct vessel wall apposition is obtained by anchoring the stent after the first rows of struts exit the microcatheter. There have been only a few reports on SAC with Neuroform Atlas stents. Most of them reported their preliminary experience with using Neuroform Atlas stents [6–13]. From their experiences with 27 consecutive patients one study concluded that Neuroform Atlas SAC is a feasible way to treat ruptured and unruptured wide-necked aneurysms that are not amenable to conventional coiling or clipping [9]. Aneurysm occlusion and favorable clinical outcomes were consistent with previously reported rates for SAC of wide-neck aneurysms using other devices [9]. From their

Park et al. BMC Neurology (2019) 19:245 Page 8 of 9

experiences with 37 aneurysms in 36 patients, another report concluded that deployment of the Neuroform Atlas stent was a safe and effective method for the treatment of wide-necked aneurysms [8].

Although it is not always easy to advance a microcatheter into an aneurysm using SAC with the cell-through technique, there are some advantages in the use of the Neuroform Atlas stent for the cell-through technique. (1) The closed cells at the proximal end with flaring end give good navigability and enhance the cell-through technique by passing a microwire in the stent lumen easily. When a microwire cannot be advanced through the stent or is stuck in the struts, a microwire looping technique (Fig. 1d and e) provides easy navigation through the inside of the stent lumen. (2) Its open-cell design with both ends flared gives good stability to the vessel wall and wall apposition. It is easy for microcatheters and microwires to pass through the inside of the stent without irritating the struts. Otherwise, microwires cannot be navigated through the inside of the stent. Microwires may go out of the stent lumen, between the struts and the vessel wall, then come back into the stent lumen again. (3) The struts are flexible and reduce stent width. When microcatheter tips cannot be advanced into the aneurysm due to being stuck in the struts, it is possible to overcome the situation by navigating a 0.014-in. microwire through the cell followed by advancement of the microcatheter over the microwire. Pushing the microcatheter very gently against the strut, withdrawing the microwire tip in the microcatheter allows the microcatheter to advance into the aneurysm while overcoming the thin and flexible struts (Fig. 1g and h). (4) It can be delivered via a 0.0165-in. microcatheter that can be easily placed intra-stent space through the pusher wire and easily navigated into an aneurysm sac due to a smaller size of caliber and a less ledge effect. Previous stents can be delivered via a 0.027- or 0.021-in. microcatheter that was difficult to be placed into an aneurysm sac due to a relatively larger size of its caliber and a ledge effect. In some situations, exchange by 0.017 or 0.0165 microcatheter was mandatory in order to use a 0.027- or 0.021-in. microcatheter. A benefit of this cell-through SAC with a single microcatheter method is the possibility to perform SAC using 5Fr-guiding or intermediated catheters, and even when using a 4Fr-intermediate catheter. Another benefit was easy navigation through very tortuous vessels. And, once the aneurysm selection was done by the cell-through technique, then the microcatheter was stuck within a cell of the stent, becoming very stable while inserting coils. Thus, patients with small-caliber, stenotic, or very tortuous parent arteries (ICA or VA) can be good candidates for this technique.

There are some other techniques to perform SAC rather than this cell-through technique. We could perform

SAC using two parent arteries with two 5Fr-guiding catheters in each parent artery. First, coil embolization with horizontal stenting could be performed crossing the Circle of Willis [14–19]. An aneurysm on the anterior communicating artery or the ICA bifurcation could be treated using a 5Fr-guding catheter each in bilateral ICAs. An aneurysm in the basilar artery aneurysm could be treated using a 5Fr-guding catheter in each ICA and VA. Second, we could perform SAC using a 5Fr-guiding catheter in bilateral VAs to treat aneurysms in the posterior circulation [20, 21]. An aneurysm on the posterior inferior cerebellar artery could be treated using bilateral VAs. However, none of these are universal techniques for SAC. Those techniques can be applied on a case by case basis. Thus, we believed that this cell-through SAC technique with a single microcatheter and a Neuroform Atlas stent might provide neurointerventionists one more option to treat wide-necked intracranial aneurysms with small-caliber, stenotic, or very tortuous parent arteries.

Other low-profile neurovascular stents have been introduced and used clinically, including the braided LVIS Jr. or LEO Baby (Balt, Montmorency, France) that are compatible with a 0.0165-in. microcatheter [22]. These devices can be placed through the same microcatheter needed for subsequent coiling. Thus, they may be used in this cell-through SAC. However, LEO Baby is not available in our country. We have experienced cell-through SAC using LVIS Jr., not intentionally with a single microcatheter but secondarily due to accidental kickback of the microcatheter while performing SAC with a jailing technique. Thus, we have not performed cell-through SAC using LVIS Jr. with the same single microcatheter.

This study was a retrospective study with a small case series so that the effects of possible selection bias cannot be excluded. However, the cases in the present study could not be done by usual SAC using two microcatheters. Another technique, such as cell-through SAC using a single microcatheter, was inevitable. Thus, it was difficult to have a true control group for evaluating the safety and efficacy of this cell-through SAC technique. In addition, a longer follow-up period is needed to evaluate the true safety and efficacy of this technique. Furthermore, it was challenging to describe the geometry and tortuosity of cerebral arteries in 3D space, because determining uniform measurement standards was difficult. In the present study, we grossly categorized the vessel conditions; parent arteries of small-caliber, with stenosis, or a very tortuous course.

# **Conclusions**

When performing SAC in patients with wide-necked intracranial aneurysms, rather than struggling to place a Park et al. BMC Neurology (2019) 19:245 Page 9 of 9

6Fr-guiding catheter in small-caliber, stenotic, or very tortuous parent arteries (ICA or VA), the cell-through SAC might be a useful option by using a single microcatheter and a Neuroform Atlas stent within a 5Fr- (or smaller) guiding or intermediate catheter.

#### **Abbreviations**

CT: Computed tomography; DSA: Digital subtraction angiography; ICA: Internal carotid artery; MRA: Magnetic resonance angiography; mRS: Modified Rankin Scale; SAC: Stent-assisted coiling; VA: Vertebral artery

#### Acknowledgements

None.

#### Authors' contributions

KYP had contributions to the concept of the study, interpretation of data, drafting the manuscript, and revising the manuscript. CKJ had contributions to the acquisition of data, analysis of data, and revising the manuscript. JWL, DJK, and BMK had contributions to the interpretation of data and revising the manuscript. JC contributed to the concept of the study, revising and editing the manuscript, and supervising the whole process of the study. All authors have approved the final version of the manuscript.

#### Funding

No funding was received for this study.

### Availability of data and materials

The datasets used and analyzed during the current study are available from the corresponding author on reasonable request.

# Ethics approval and consent to participate

All procedures performed in the studies involving human participants were in accordance with the ethical standards of our Institutional Review Board with the 1964 Helsinki Declaration and its later amendments or comparable ethical standards. In this retrospective study, the requirement for informed consent was waived.

### Consent for publication

Not applicable.

# Competing interests

The authors declare that they have no competing interests.

#### **Author details**

<sup>1</sup>Department of Neurosurgery, Severance Hospital, Yonsei University College of Medicine, 50-1, Yonsei-ro, Seodaemun-gu, Seoul 03722, Republic of Korea. <sup>2</sup>Department of Radiology, Severance Hospital, Yonsei University College of Medicine, Seoul, Republic of Korea. <sup>3</sup>Severance Institute for Vascular and Metabolic Research, Yonsei University College of Medicine, Seoul, Republic of Korea.

# Received: 2 July 2019 Accepted: 20 September 2019 Published online: 22 October 2019

# References

- Pierot L, Biondi A. Endovascular techniques for the management of wideneck intracranial bifurcation aneurysms: a critical review of the literature. J Neuroradiol. 2016;43(3):167–75.
- Chung J, Lim YC, Suh SH, et al. Stent-assisted coil embolization of ruptured wide-necked aneurysms in the acute period: incidence of and risk factors for periprocedural complications. J Neurosurg. 2014;121(1):4–11.
- Tahtinen Ol, Vanninen RL, Manninen HI, et al. Wide-necked intracranial aneurysms: treatment with stent-assisted coil embolization during acute (
  72 hours) subarachnoid hemorrhage-experience in 61 consecutive patietns. Radiology. 2009;253(1):199–208.
- Mascitelli JR, Moyle H, Oermann EK, et al. An update to the Raymond–Roy occlusion classification of intracranial aneurysms treated with coil embolization. J Neurointerv Surg. 2015;7(7):496–502.

- Park MS, Stiefel MF, Fiorella D, Kelly M, McDougall CG, Albuquerque FC. Intracranial placement of a new, compliant guide catheter: technical note. Neurosurgery. 2008;63(3):E616–7.
- Putman CM, Chaloupka JC. Use of large-caliber coronary guiding catheters for neurointerventional applications. AJNR Am J Neuroradiol. 1996;17(4): 697–704
- Turk A, Manzoor MU, Nyberg EM, Turner RD, Chaudry I. Initial experience with distal guide catheter placement in the treatment of cerebrovascular disease: clinical safety and efficacy. J Neurointerv Surg. 2013;5(3):247–52.
- Ulfert C, Pham M, Sonnberger M, et al. The Neuroform atlas stent to assist coil embolization of intracranial aneurysms: a multicentre experience. J Neurointerv Surg. 2018;10(12):1192–6.
- Ten Brinck MFM, de Vries J, Bartels RHMA, Grotenhuis JA, Boogaarts HD. Neuroform atlas stent-assisted coiling: preliminary results. Neurosurgery. 2019;84(1):179–89.
- Jankowitz BT, Hanel R, Jadhav AP, et al. Neuroform atlas stent system for the treatment of intracranial aneurysm: primary results of the atlas humanitarian device exemption cohort. J Neurointerv Surg. 2019;11(8):801–6.
- 11. Cay F, Peker A, Arat A. Stent-assisted coiling of cerebral aneurysms with the Neuroform atlas stent. Interv Neuroradiol. 2018;24(3):263–9.
- Goertz L, Dorn F, Siebert E, et al. Safety and efficacy of the Neuroform atlas for stent-assisted coiling of intracranial aneurysms: a multicenter experience. J Clin Neurosci. 2019[Epub ahead of print]. https://doi.org/10.1016/j.jocn. 2019.07.030.
- Gross BA, Ares WJ, Ducruet AF, Jadhav AP, Jovin TG, Jankowitz BT. A clinical comparison of atlas and LVIS Jr stent-assisted aneurysm coiling. J Neurointerv Surg. 2019;11(2):171–4.
- Ahmed ME, Lum C, Lesiuk H, Lancu D, dos Santos M. Navigation of stents across communicating arteries for aneurysm embolization. Can J Neurol Sci. 2014;41(2):193–9.
- Blackburn SL, Kadkhodayan Y, Shekhtman E, Derdeyn CP, Cross DT 3rd, Moran CJ. Treatment of basilar tip aneurysms with horizontal PCA to PCA stent-assisted coiling: case series. J Neurointerv Surg. 2013;5(3):212–6.
- Cho YD, Kim KM, Lee WJ, Kang HS, Kim JE, Han MH. Retrograde stenting through the posterior communicating artery in coil embolization of the posterior communicating artery aneurysm. Neuroradiology. 2013;55(6):733–9.
- Kim YJ, Sung JH, Hong JT, Lee SW. Horizontal stent assisted coiling of wide neck basilar tip aneurysm: comparison of two cases. J Cerebrovasc Endovasc Neurosurg. 2017;19(3):201–6.
- Saatci I, Geyik S, Yavuz K, Cekirge S. X-configured stent-assisted coiling in the endovascular treatment of complex anterior communicating artery aneurysms: a novel reconstructive technique. AJNR Am J Neuroradiol. 2011; 32(6):E113–7.
- Siddiqui MA, Bhattacharya JJ, Lindsay KW, Jenkins S. Horizontal stentassisted coil embolization of wide-necked intracranial aneurysms with the Enterprise stent – a case series with early angiographic follow-up. Neuroradiology. 2009;51(6):411–8.
- Cho YD, Kang HS, Lee WJ, Kim KM, Kim JE, Han MH. Stent-assisted coil embolization of wide-necked posterior inferior cerebellar artery aneurysms. Neuroradiology. 2013;55(7):877–82.
- Kim MJ, Chung J, Kim SL, et al. Stenting from the vertebral artery to the posterior inferior cerebellar artery. AJNR Am J Neuroradiol. 2012;33(2):348–52.
- Aydin K, Arat A, Sencer S, Barburoglu M, Men S. Stent-assisted coiling of wide-necked intracranial aneurysms using low-profile LEO baby stents: initial and midterm results. AJNR Am J Neuroradiol. 2015;36(10):1934–41.

# **Publisher's Note**

Springer Nature remains neutral with regard to jurisdictional claims in published maps and institutional affiliations.